

### **Protocol for AAV production**

#### **General Remarks:**

This protocol is designed for the production of two AAViruses in parallel. You need in total three 10-stacks and three 1-stacks for two rounds of transfection and harvesting. Note that two 10/1-stacks are re-used in a second round of production of the same virus (mark the stacks). The third 10/1-stack is for propagation of 293T cells.

#### **Step 1: Cell-stack expansion**

Materials required:

6 confluent (~90%) 15cm dishes 293T cells  
10-stack, 1-stack  
Trypsin (8ml)  
Medium w/o FCS 100ml or PBS  
Medium 1200ml

Cells from 6 confluent (~90%) 15cm dishes are plated onto 1x 10-stack and 1x 1-stack (this corresponds to a ~1:10 splitting).

- Following trypsinization and centrifugation cells are distributed equally to 2x 550ml medium (2 medium bottles + 10%FCS + 1%P/S).

- 2x 500ml are transferred to a 10-stack; 2x 50ml are transferred to a 1-stack.

(This is generally how a cell suspension from 2 medium bottles are distributed to a 10-stack and the corresponding 1-stack (10/1-stacks) throughout the whole protocol in order to ensure, that the cell density in the 1-stack represents the 10-stack)

## **Step 2: Pre-transfection**

### Day -1/-2, a.m.

Materials required for 2x viruses:

Cells: 1x 10-stack, 1x 1-stack 80-100% confl.  
2x 10-stacks, 2x 1-stacks  
2x 1lt Erlenmeyer flasks  
2x conical 500ml tube  
220ml 10x D-PBS (Invitrogen)  
2lt Braun- H<sub>2</sub>O (DKFZ)  
120ml Trypsin(0.25%)/EDTA (Invitrogen)  
3.3lt medium (6 bottles+FCS, P/S)

### For 2 viruses:

- Remove medium from 10-/1-stack by decanting into an Erlenmeyer flask
- Wash with 900ml/100ml D-PBS, decant into flask (minimize detachment of cells by indirectly introducing of D-PBS hitting the wall of the stack first).
- Remove D-PBS; add 100ml/10ml trypsin(0.25%)/EDTA, 37°C, 5min
- Swirl and slam the Cell-Stacks to detach cells.
- Add 500ml/50ml medium (incl. FCS) and transfer to 2x500ml conical tubes.
- Spin 500g (2000rpm), 10min, 15°C
- In the meanwhile: wash "old/used" 10-/1-stack with 900ml/100ml D-PBS
- Resuspend cells from both conical tubes into 110ml medium. For this, use 25ml from each of 4 medium bottles (for 2 10/1-stacks to be transfected for the production of 2 viruses); and 5ml from 2 medium bottles (for the third, already used and washed 10/1-stack for propagation of cells for the second round of transfection).
- Distribute the cell suspension to the according medium bottles:  
4x25ml (transfection); 2x5ml (propagation)
- Fill the three 10/1-stacks with 1000ml/100ml of the corresponding cell suspensions.

### Step 3: Transfection

#### Day 0, a.m.

Materials required for 2x viruses ( 2x10-stacks+2x1-stacks):

395µg dsAAV vector (for each virus construct)

994µg rep/cap plasmid (p5E18-VD2/8)

2706µg pDGΔVP

40ml H<sub>2</sub>O (UltraPure)

36ml PEI (7,5 mM; see protocol for production of PEI)

80ml NaCl (300mM)

4x sterile bottles (min. 100ml; for mixing of transfection reagent and DNA)

2x conical 500ml tubes

- Cells should be 70-80% confluent

- Prepare DNA solution: **Mix1**

395µg pdsAAV vector (~5.8kb)

497µg p5E18-VD2/8 (~7.3kb) (for AAV2/8)

1353µg pDGΔVP (~21kb)

H<sub>2</sub>O ad 20 ml

+ 20ml NaCl (300mM); mix

40ml

- Prepare PEI solution: **Mix2**

18ml PEI

2 ml H<sub>2</sub>O

+ 20ml NaCl (300mM); mix

40ml

- Introduce mix 2 dropwise into mix 1 while swirling (immediate mixing avoids the generation of large DNA precipitates that are not transfectable; solution gets cloudy)
- Incubate for 10 minutes at room temperature
- Transfer 40ml of transfection mix to each of the two medium bottles/virus (including FCS and P/S)
- Remove "old" medium from the 1- and the 10-stack
- Transfer 50ml from each medium bottle to the 1-stack (100ml medium) and the remaining medium to the 10-stack (containing ~2x550ml medium).
- Change medium 24h post-transfection.

#### **Step 4: Harvest**

Day 2, p.m. (or day 3, a.m.)

Materials required for 2x viruses:

> 2 1lt-Erlenmeyer flasks  
550ml 10x D-PBS  
5lt Braun-H<sub>2</sub>O  
230ml Trypsin(0.25%)/EDTA  
800ml medium (incl. FCS)  
4x conical 500ml tubes  
Lysis buffer (150mM NaCl, 50mM TrisHCl pH8.5, sterile)

For each 10/1-Stack (virus):

- Remove medium by decanting into Erlenmeyer flasks
- Wash with 900ml/100ml D-PBS, decant into flask (if necessary spin detached cells down @ 600g, 10min, 15°C)
- Add 100ml/10ml trypsin(0.25%)/EDTA, 37°C, 5min
- Swirl and slam the Cell-stack to detach cells
- Add 350ml/40ml medium (incl. FCS) and transfer to a 500ml conical tube (390ml + 110 ml trypsin = 500 ml); (1<sup>st</sup> tube).
- Add 450ml/50ml D-PBS and transfer to a 500ml conical tube (2<sup>nd</sup> tube).
- Spin both tubes 600g, 10min, 15°C
- In the meantime Wash Cell-stacks with 500ml/50ml D-PBS (1<sup>st</sup> round only) and add 500ml/50ml D-PBS for storage (1-2 days in 4°C room for 2<sup>nd</sup> round).
- 2<sup>nd</sup> tube: Discard SN, resuspend pellet in 50ml D-PBS and transfer to 50ml-Falcon
- 1<sup>st</sup> tube: Transfer SN to 2<sup>nd</sup> conical tube, use 50ml cell suspension from 2<sup>nd</sup> Falcon to resuspend pellet and distribute to 2x 50ml-Falcons.
- Spin Falcons 2000rpm, 10min, 4°C (swing-rotor centrifuge)
- Resuspend pellets in 8ml lysis buffer and transfer one 50ml-Falcon
- Vortex and snap-freeze in N<sub>2</sub>, store at -80°C.

---

Repeat Step 2 to 4 for a second round of virus production reusing the 10/1-stacks for the same virus as before.

### Step 5: Lysis (cont.)

Materials required for 2x viruses:

N<sub>2</sub>

Lysis buffer

Sterile MgCl<sub>2</sub> 1M

Benzonase >8000U (SIGMA)

Water bath sonicator

Starting with 2x50ml Falcons per virus containing 2x~20ml lysate from two rounds of transfection:

- Thaw in water bath at 37°C, vortex and snap-freeze in N<sub>2</sub>
- Repeat once
- Thaw, vortex, transfer into fresh 4x15ml-Falcons and spin ~3500g (6000 rpm TC centrifuge), 10min, 4°C
- Collect **SN1** (4x~10ml) into 4x15ml-Falcons and resuspend pellet in 5ml lysis buffer
- Vortex, snap-freeze, thaw at 37°C
- Repeat twice
- Transfer into fresh 4x15ml-Falcon and spin ~3500g, 10min, 4°C
- Collect **SN2** (4x~10ml) into 4x15ml-Falcons and resuspend pellet in 4ml lysis buffer
- Vortex, snap-freeze, thaw at 37°C.
- Vortex, sonicate in water bath sonicator 1min.
- Add MgCl<sub>2</sub> to 5mM and Benzonase (Sigma) to 50U/ml to **SN1+2** and to lysate (3)
- Incubate at 37°C for 30min
- Transfer lysate (3) to fresh 15ml-Falcon and spin ~3500g, 10min, 4°C
- Collect **SN3** and pool with **SN2**
- Freeze 4xSN1 (~10ml) and 4xSN2/3 (~9ml) and 4xpellets in (-80°C)

Optional: take 20µl aliquots from each fraction for possible qPCR tests (for troubleshooting in case virus yield is low at the end.)

## **Step 6: Gradient purification**

### Materials required for 2 viruses:

10x Beckman Centrifuge Tubes, Polyallomer, Quick-Seal, 39ml, 25x89  
4x spacers for centrifuge tubes  
UZ, 50.2Ti rotor, 3x 3 hours  
126ml OptiPrep (60% iodixanol)  
270µl phenol red 0.5%  
MgCl<sub>2</sub> 1M, sterile  
KCl 4M, sterile  
NaCl 5M, sterile  
120ml D-PBS 1x (Invitrogen)  
Long pasteur pipettes  
Lysis buffer  
50ml syringes  
10ml syringes  
5ml syringes  
filters for 50ml syringes  
20G needles

### For 1 virus (4 gradients):

- Prepare iodixanol dilutions while thawing lysates.

For 2 viruses prepare the according amounts (see next page) for 10 gradients for the 1<sup>st</sup> round of centrifugation and 5 gradients for the 2<sup>nd</sup> round.

- Pool each SN1 with SN2/3: 4x~19ml:
- Transfer lysates to tubes through long Pasteur pipette.
- Sequentially underlayer through Pasteur pipette (use prefilled syringes for all layers of one gradient to ensure continuous (air-bubble free) filling of the tube):

7ml 15% iodixanol dil.

5ml 25% iodix. dil.

4ml 40% iodix. dil.

4ml 60% iodix. dil.

- Top up and balance pair-wise with lysis buffer so that ~1ml air remains on top
- Seal tube nozzle with heat and place spacers
- Load ultracentrifuge in 50.2Ti rotor and set up vacuum
- Spin at 50krpm, 2.5h, 10°C
- Insert a 20G needle on top of tube and into the lysate layer

## AAV Cellstack production

A. Vegiopoulos; modified: M.Berriel Diaz 12/2011

- Insert a 5ml syringe with 20G needle into the 60% layer and collect max. 3.5ml from the 40% layer from underneath. Avoid collecting from the 25% layer
- Separate aliquots as required (20µl), pool suspensions from 4 gradients per virus (~14ml; can be stored at -80°C).

Iodixanol dilution	1 gradient	3 gradients	5 gradients	10 gradients	12 gradients
15%	7ml	21ml	35ml	70ml	84ml
25%	5ml	15ml	25ml	50ml	60ml
40%	4ml	12ml	20ml	40ml	48ml
60%	4ml	12ml	20ml	40ml	48ml

15% iodix. dil.	1 gradient	3 gradients	5 gradients	10 gradients	12 gradients
OptiPrep	1.75ml	5.25ml	8.75ml	17.5ml	21ml
PBS-MK-NaCl	5.25ml	15.75ml	26.25ml	52.5ml	63ml
Final	7ml	21ml	35ml	70ml	84ml

25% iodix. dil.	1 gradient	3 gradients	5 gradients	10 gradients	12 gradients
OptiPrep	2.08ml	6.24ml	10.4ml	20.8ml	25ml
PBS-MK	2.92ml	8.76ml	14.6ml	29.2ml	35ml
Phenol red 0.5%	12.5µl	37.5µl	62.5µl	125µl	150µl
Final	5ml	15ml	25ml	50ml	60ml

40% iodix. dil.	1 gradient	3 gradients	5 gradients	10 gradients	12 gradients
OptiPrep	2.67ml	8ml	13.35ml	26.7ml	32ml
PBS-MK	1.33ml	4ml	6.65ml	13.3ml	16ml
Final	4ml	12ml	20ml	40ml	48ml

60% iodix. dil.	1 gradient	3 gradients	5 gradients	10 gradients	12 gradients
OptiPrep	4ml	12ml	20ml	40ml	48ml
Phenol red 0.5%	10µl	30µl	50µl	100µl	120µl
Final	4ml	12ml	20ml	40ml	48ml

PBS-MK*	Final	1 gradient	3 gradients	5 gradients	10 gradients	12 gradients
D-PBS 1x	1x	6ml	18ml	24ml	45ml	55ml
MgCl <sub>2</sub> 1M	1mM	6µl	18µl	24µl	45µl	55µl
KCl 4M	2.5mM	3.75µl		15µl	28µl	34.37µl

\* sterilize by filtration

PBS-MK-NaCl*	Final	1 grad.	3 grad.	5 grad.	10 grad.	12 grad.
D-PBS 1x	~0.8x	5.6ml	16.8ml	23.2ml	44ml	52.8ml
MgCl <sub>2</sub> 1M	1mM	7µl	21µl	29µl	55µl	66µl
KCl 4M	2.5mM	4.37µl	13.11µl	18.12µl	34µl	41.25µl
NaCl 5M	1M	1.4ml	4.2ml	5.8ml	11ml	13.2ml

\* sterilize by filtration

### 2<sup>nd</sup> round gradients

- Use 1-2 gradients per virus
- Add lysis buffer to suspensions up to 25ml total volume and transfer to a Beckman tube
- Sequentially underlayer with:
  - 5ml 25% iodix. dil.
  - 4ml 40% iodix. dil.
  - 4ml 60% iodix. dil.
- Top up and balance pair-wise with lysis buffer so that ~1ml air remains on top
- Seal tube nozzle with heat (Kleinschmidt) and place spacers
- Load ultracentrifuge in 50.2Ti rotor and set up vacuum
- Spin at 50krpm, 2.5h, 10°C
- Insert a 20G needle on top of tube and into the lysate layer
- Insert a 5ml syringe with 20G needle into the 60% layer and collect max. 3.5ml from the 40% layer from underneath. Avoid collecting from the 25% layer
- Separate aliquots and proceed to dialysis or freeze at -80°C

### **Dialysis and concentration**

#### Materials required for 2 viruses:

Slide-A-Lyzer-kit 10 KD 3-12 mL (Perbio)  
 300ml 10x D-PBS  
 3ltr Braun-H<sub>2</sub>O

AAV Cellstack production

A. Vegiopoulos; modified: M.Berriel Diaz 12/2011

- Dialyse in slide a lycer (12ml extra strength) following instructions in the kit:
- Dialyse against 1l D-PBS each for  
3h at 4°C  
over night at 4°C  
3h at 4°C
- Transfer into 15 ml Falcon tubes.
- Transfer gradually into Vivaspin tubes for concentration; spin at 2000-4000rpm for 3-5 min until the final volume is about 500-1000µl; piped up and down while loading new volume on the column to remove virus from the filter walls.
- Take a 5µl aliquot for titration.
- Make 200-250µl aliquots and freeze at -80°C.