

Monochromatic Imaging in Dual Energy CT (DECT): Metal Artifact Reduction with Acceptable Image Quality?

Stefan Kuchenbecker¹, Sebastian Faby¹, Sören Schüller¹,
Michael Lell², and Marc Kachelrieß^{1,3}




¹German Cancer Research Center (DKFZ), Heidelberg, Germany

²Universitätsklinikum, Erlangen, Germany

³Friedrich-Alexander-University (FAU), Erlangen, Germany

DECT and Pseudo Monochromatic Imaging

Pseudo monochromatic imaging is a linear combination of DECT f_L and f_H : $f_\alpha = (1 - \alpha) f_L + \alpha f_H$



Monochromatic Imaging

- **Pseudo monochromatic imaging** $f_\alpha = (1 - \alpha) f_L + \alpha f_H$
 - Image-based postprocessing
 - Provided in clinical DECT scanners
- **Virtual monochromatic imaging** $g_\alpha = (1 - \alpha) g_L + \alpha g_H$
 - Rawdata-based preprocessing
 - Constraint on consistent rawdata
- **True monochromatic imaging**
 - Would require monochromatic x-rays – not applicable here


$$q_L = -\ln \int dE w_L(E) e^{-p_W \mu_W(E) - p_B \mu_B(E)}$$

$$q_H = -\ln \int dE w_H(E) e^{-p_W \mu_W(E) - p_B \mu_B(E)}$$

Series Expansion

- Series expansion of the polychromatic attenuation:



$$q_j = -\ln \int dE w_j(E) e^{-p_W \mu_W(E) - p_B \mu_B(E)} = \sum_{kl} c_{jkl} p_W^k p_B^l$$



pseudo monochromatic
image-based processing



virtual monochromatic
rawdata-based processing

$f_L = f_0$
($E = 67$ keV)





g_0
($E = 67$ keV)

maximum CNR





maximum CNR

$f_H = f_1$
($E = 93$ keV)



g_1
($E = 93$ keV)

$f_{1.67}$
($E = 221$ keV)



$g_{1.67}$
($E = 221$ keV)



(40 / 400)

Patient Data Set – Pseudo Monochromatic Imaging


$z = -178$ mm

$z = -264$ mm


$f_L = f_0$
($E = 67$ keV)





$f_H = f_1$
($E = 93$ keV)



$f_{1.50}$
($E = 140$ keV)






$f_{1.90}$
($E = \dots$ keV)



$C = 0$ HU, $W = 800$ HU

Conclusion

- Pseudo monochromatic imaging
 - cannot completely remove metal artifacts,
 - can sometimes reduce metal artifacts,
 - reduces CNR, if used for metal artifact reduction.
- Rawdata-based DECT decomposition is to be preferred.
- Rawdata-based MAR methods such as FSNMAR^{1,2} should be preferred.
- The additional information available in DECT should be used for spectral imaging rather than for artifact reduction.



Thank You!

This study was supported by the Deutsche Forschungsgemeinschaft under grants KA 1678/5-1 and LE 2763/1-1.

Parts of the reconstruction software were provided by RayConStruct® GmbH, Nürnberg, Germany.

This presentation will soon be available at www.dkfz.de/ct.